BladderScan i10™

BladderScan® Prime Plus™

- EN** Calibration Manual
- DA** Kalibreringsmanual
- DE** Kalibrierungshandbuch
- ES** Manual de calibración
- FI** Kalibrintiopus

- FR-CA** Manuel d'étalonnage
- FR-FR** Manuel d'étalonnage
- IT** Manuale di calibrazione
- NL** Kalibratiehandleiding
- NB** Kalibreringshåndbok
- SV** Kalibreringsmanual

BladderScan
verathon

BladderScan i10™

BladderScan® Prime Plus™

- EN** Calibration Manual
- DA** Kalibreringsmanual
- DE** Kalibrierungshandbuch
- ES** Manual de calibración
- FI** Kalibrointiopas
- FR-CA** Manuel d'étalonnage
- FR-FR** Manuel d'étalonnage
- IT** Manuale di calibrazione
- NL** Kalibratiehandleiding
- NB** Kalibreringshåndbok
- SV** Kalibreringsmanual

Effective: September 21, 2022

Caution: Federal (United States) law restricts this device to sale by or on the order of a physician.

Contact Information

To obtain additional information regarding your BladderScan system, please contact Verathon Customer Care or visit [verathon.com/service-and-support](https://www.verathon.com/service-and-support).

Verathon Inc.

20001 North Creek Parkway
Bothell, WA 98011 U.S.A.
Tel: +1 800 331 2313 (US and Canada only)
Tel: +1 425 867 1348
Fax: +1 425 883 2896
[verathon.com](https://www.verathon.com)

Verathon Medical (Europe) B.V.

Willem Fenengastraat 13
1096 BL Amsterdam
The Netherlands
Tel: +31 (0) 20 210 30 91
Fax : +31 (0) 20 210 30 92

Verathon Medical (Australia) Pty Limited

Unit 9, 39 Herbert Street
St Leonards NSW 2065
Australia
Within Australia: 1800 613 603 Tel / 1800 657 970 Fax
International: +61 2 9431 2000 Tel /
+61 2 9475 1201 Fax

MDSS CH GmbH

Laurenzenvorstadt 61
5000 Aarau
Switzerland

CH

Anandic Medical Systems AG

Stadtweg 24
8245 Feuerthalen
Switzerland

MDSS-UK RP Ltd.

6 Wilmslow Road, Rusholme
Manchester M14 5TP
United Kingdom
Tel: +44 (0)7898 375115

Copyright © 2022 Verathon, Inc. All rights reserved.

BladderScan, BladderScan i10, BladderScan Prime Plus, the BladderScan symbol, and Verathon are trademarks of Verathon Inc. All other brand and product names are trademarks or registered trademarks of their respective owners.

Information in this manual may change at any time without notice. For the most up-to-date information, see the documentation available at [verathon.com/service-and-support](https://www.verathon.com/service-and-support).

California residents — For Proposition 65 information, visit [verathon.com/Prop65](https://www.verathon.com/Prop65).

Table of Contents

- EN** English.....1
- DA** Dansk.....7
- DE** Deutsch.....13
- ES** Español.....19
- FI** Suomi25
- FR-CA** Français (Canada).....31
- FR-FR** Français (France)37
- IT** Italiano43
- NL** Nederlands47
- NB** Norsk.....53
- SV** Svenska59


Calibrating the System

If your facility requires it, you can calibrate your BladderScan system periodically as described in the following procedure. The BladderScan system does not require routine or periodic calibration under normal use. For all other information about how to use and maintain your BladderScan system, see its Operations and Maintenance Manual, which is available at verathon.com/service-and-support.

Procedure 1. Perform the Initial Inspection

To calibrate the BladderScan system, you must have a calibration tank designed for it. When you receive the calibration tank, be sure to inspect it before you use it.

1. Carefully open the top of the shipping box. Do not insert anything sharp through the box.
2. Remove the contents and verify that you have received the components in the following figure.


3. Inspect all components for damage.
4. If any of the components are missing or damaged, notify the carrier and Verathon Customer Care or your local representative.


Procedure 2. Perform the Calibration (BladderScan i10 Systems)

If a PIN for Settings has been defined in the BladderScan i10 system, you will need to know that PIN to complete this procedure. However, no unique authorization PIN is required.

1. Place the calibration tank on a flat, non-reflective surface, and then remove the lid.
2. Install the spiral-shaped target in the tank base, using the notches in the base to position the target correctly.
3. Pour clean, room-temperature water into the container, filling to the top of the textured section. Add the water slowly, without splashing, to avoid mixing air into it.
4. Allow the water in the calibration tank to sit for at least 24 hours so that any embedded air can escape.
5. After the air has settled out of the water, skim or tap the target and the tank to remove as many air bubbles as possible.


6. Replace the tank lid onto the base.
7. Place the probe into the collar in the tank lid. Align the probe button with the notch exactly as shown, and then verify that the tip of the probe is submerged in the water.

8. If necessary, tap or agitate the probe to remove as many air bubbles as possible from its submerged tip.
9. On the system's Home screen, tap **Settings**, and then tap **Admin**.
10. If a keypad appears, key in the PIN for Settings, and then tap **Enter**.

11. At the bottom of the screen, tap **Calibrate**.

12. Review the instructions that appear on the screen, and then tap the **Calibrate** button above them. A waiting indicator appears on the screen while the calibration is in progress.

*Note: Once the calibration process has started, do not press the **On/Off** button or disconnect the probe from the console.*

- If calibration is successful, a message on the console indicates that the probe has passed, and the Last Calibration and Calibration Expiration dates are updated.
- If calibration fails, a message on the console indicates the failure.

13. If a calibration failure occurs:

- Make sure the calibration tank contains enough water.
- Remove as many air bubbles as possible from the water, the target, and the probe.
- Verify that the probe is seated properly in the tank lid.

14. If you want to calibrate the system again, tap **Calibrate** to begin recalibration. If calibration fails repeatedly, contact Verathon Customer Care.


15. When calibration is complete, tap **Back**.

16. Remove the probe from the tank lid, and then dry it with a clean, soft cloth.


Procedure 3. Perform the Calibration (BladderScan Prime Plus Systems)

You do not need to obtain a unique authorization PIN to calibrate a BladderScan Prime Plus system.

1. Place the calibration tank on a flat, non-reflective surface, and then remove the lid.
2. Install the spiral-shaped target in the tank base, using the notches in the base to position the target correctly.
3. Pour clean, room-temperature water into the container, filling to the top of the textured section. Add the water slowly, without splashing, to avoid mixing air into it.
4. Allow the water in the calibration tank to sit for at least 24 hours so that any embedded air can escape.
5. After the air has settled out of the water, skim or tap the target and the tank to remove as many air bubbles as possible.


6. Replace the tank lid onto the base.
7. Place the probe into the collar in the tank lid. Align the probe button with the notch exactly as shown, and then verify that the tip of the probe is submerged in the water.

8. If necessary, tap or agitate the probe to remove as many air bubbles as possible from its submerged tip.
9. On the system's Home screen, tap **Settings**.

10. Tap **Admin**, and then at the bottom of the screen, tap **Calibrate**.

11. Review the instructions that appear on the screen, and then tap the **Calibrate** button above them. A waiting indicator appears on the screen while the calibration is in progress.

*Note: Once the calibration process has started, do not press the **On/Off** button or disconnect the probe from the console.*

- If calibration is successful, a message on the console indicates that the probe has passed, and the Last Calibration and Calibration Expiration dates are updated.
- If calibration fails, a message on the console indicates the failure.

12. If a calibration failure occurs:

- Make sure the calibration tank contains enough water.
- Remove as many air bubbles as possible from the water, the target, and the probe.
- Verify that the probe is seated properly in the tank lid.

13. If you want to calibrate the system again, tap **Calibrate** to begin recalibration. If calibration fails repeatedly, contact Verathon Customer Care.

14. When calibration is complete, tap **Back**.

15. Remove the probe from the tank lid, and then dry it with a clean, soft cloth.

Procedure 4. Disassemble the Calibration Tank (Optional)

Follow these steps if you store the calibration tank between uses.

1. Remove the tank lid from the tank base.
2. Remove the target from the base.
3. Empty the water from the base.
4. Using a soft, clean cloth, dry the target, the base, and the lid.

CAUTION

If you empty the calibration tank and store it between uses, store it in a location that protects it from impact and other damaging external forces.

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