


# Single-Use Without a Single Compromise



## GlideScope Spectrum QC eco

New generation Spectrum QC eco™ video laryngoscopes couple innovative magnetic QuickConnect™ technology and imaging improvements with new advances towards sustainability. These enhancements take Spectrum™ video laryngoscopes to a new level of capability—for patients, physicians and the planet.



### Innovative Connections


Magnetic QuickConnect for plug-and-play simplicity

### Fast, Efficient Care

QuickConnect for Spectrum QC and BFlex™ simplifies your GlideScope® Core™ system for procedural efficiency



Broad portfolio of nine blade sizes and types



### Enhanced Airway Views

50% improvement in airway image sharpness<sup>1</sup>









### Expanded Visualization

14% greater field of view expands visualization in the airway<sup>2</sup>

1. Comparing Spectrum QC Hyperangle/Mac 3 & 4 versus Spectrum HDMI LoPro/Mac 3 & 4

2. Diagonal field of view, comparing Spectrum QC Hyperangle/Mac versus Spectrum HDMI LoPro/Mac

# Spectrum Product Comparison

Workflow and Connectivity	Spectrum QC eco™	Spectrum™
Connector Type	Magnetic QuickConnect™ 	HDMI-like 
Connector Retention Force	50% higher retention force <sup>1</sup>	–
Dual-View Cable Requirements	Two GlideScope Core QuickConnect Cables	One Core Smart Cable One Core QuickConnect Cable
<b>Imaging</b>		
Field of View	14% greater field of view <sup>2</sup>	–
Dynamic Light Control®	Yes	Yes
LED Brightness	Mac 30% brighter; Miller 100% brighter <sup>3</sup>	–
Image Sharpness	50% sharper <sup>4</sup>	–
<b>Visualization System Compatibility</b>		
GlideScope® Core™ Monitors 	GlideScope Core QuickConnect Cable (for Spectrum QC and BFlex™) 	GlideScope Core Smart Cable (for Spectrum™) 
GlideScope® Go™ 	GlideScope Go 2	GlideScope Go
<b>Sustainability Advancements</b>		
Bio-Based Plastic Shell <sup>5</sup>	Yes	No
Avoids EtO Gas Sterilization <sup>6</sup>	Yes	No
Recycled Shipment Packaging	Yes	Yes

1. Spectrum QC with Core QuickConnect cable versus Spectrum with Core Smart Cable
2. Diagonal field of view, comparing Spectrum QC Hyperangle/Mac versus Spectrum HDMI LoPro/Mac
3. Comparing Spectrum QC Mac/Miller versus Spectrum Mac/Miller
4. Spectrum QC Hyperangle/Mac 3 & 4 versus Spectrum LoPro/Mac 3 & 4
5. The raw material of the shells is linked to 80% bio-circular feedstock, which is allocated via the mass balance approach based on the amount of bio-circular material sourced in its production.
6. Spectrum QC eco is non-sterile. It is produced and packaged in an ISO Class 8 cleanroom environment, ready for patient use.

To learn more, contact your local Verathon representative or visit [verathon.com/spectrumqc](https://verathon.com/spectrumqc)



[verathon.com](https://verathon.com)

Verathon Inc.  
20001 North Creek Parkway  
Bothell, WA 98011  
USA  
Tel: +1 800 331 2313  
(USA & Canada only)  
Tel: +1 425 867 1348  
Fax: +1 425 883 2896



Verathon Medical (Europe) B.V.  
Willem Fenengastrat 13  
1096 BL Amsterdam  
The Netherlands  
Tel: +31 (0) 20 210 30 91  
Fax: +31 (0) 20 210 30 92

Verathon Medical (Australia) Pty Limited  
Unit 9, 39 Herbert Street  
St Leonards NSW 2065  
Australia  
Tel: 1800 613 603 (Within Australia)  
Fax: 1800 657 970  
Tel: +61 2 9431 2000 (International)  
Fax: +61 2 9475 1201