


Single-Use Without a Single Compromise



GlideScope Spectrum QC eco

New generation Spectrum QC eco™ video laryngoscopes couple innovative magnetic QuickConnect™ technology and imaging improvements with new advances towards sustainability. These enhancements take Spectrum™ video laryngoscopes to a new level of capability—for patients, physicians and the planet.




Innovative Connections


Magnetic QuickConnect for plug-and-play simplicity

Fast, Efficient Care

QuickConnect for Spectrum QC and BFlex™ simplifies your GlideScope® Core™ system for procedural efficiency



Broad portfolio of thirteen blade sizes and types



Enhanced Airway Views

50% improvement in airway image sharpness¹









Expanded Visualization

14% greater field of view expands visualization in the airway²

1. Comparing Spectrum QC Hyperangle/Mac 3 & 4 versus Spectrum HDMI LoPro/Mac 3 & 4

2. Diagonal field of view, comparing Spectrum QC Hyperangle/Mac versus Spectrum HDMI LoPro/Mac

Spectrum Product Comparison

| Workflow and Connectivity | Spectrum QC eco™ | Spectrum™ |
|--|---|---|
| Connector Type | Magnetic QuickConnect™  | HDMI-like  |
| Connector Retention Force | 50% higher retention force ¹ | – |
| Dual-View Cable Requirements | Two GlideScope Core QuickConnect Cables | One Core Smart Cable and One Core QuickConnect Cable |
| Blade Options | 13 blade sizes/types | 9 blade sizes/types |
| Imaging | | |
| Field of View | 14% greater field of view ² | – |
| Dynamic Light Control® | Yes | Yes |
| LED Brightness | Mac 30% brighter; Miller 100% brighter ³ | – |
| Image Sharpness | 50% sharper ⁴ | – |
| Visualization System Compatibility | | |
| GlideScope® Core™ Monitors  | GlideScope Core QuickConnect Cable (for Spectrum QC and BFlex™)  | GlideScope Core Smart Cable (for Spectrum™)  |
| GlideScope® Go™  | GlideScope Go 2 | GlideScope Go |
| Sustainability Advancements | | |
| Bio-Based Plastic Shell ⁵ | Yes | No |
| Avoids EtO Gas Sterilization ⁶ | Yes | No |
| Recycled Shipment Packaging | Yes | Yes |

1. Spectrum QC with Core QuickConnect cable versus Spectrum with Core Smart Cable
2. Diagonal field of view, comparing Spectrum QC Hyperangle/Mac versus Spectrum HDMI LoPro/Mac
3. Comparing Spectrum QC Mac/Miller versus Spectrum Mac/Miller
4. Spectrum QC Hyperangle/Mac 3 & 4 versus Spectrum LoPro/Mac 3 & 4
5. The raw material of the shells is linked to 80% bio-circular feedstock, which is allocated via the mass balance approach based on the amount of bio-circular material sourced in its production.
6. Spectrum QC eco is non-sterile. It is produced and packaged in an ISO Class 8 cleanroom environment, ready for patient use.

To learn more, contact your local Verathon representative or visit verathon.com/spectrumqc



verathon.com

Verathon Inc.
20001 North Creek Parkway
Bothell, WA 98011
USA
Tel: +1 800 331 2313
(USA & Canada only)
Tel: +1 425 867 1348
Fax: +1 425 883 2896



Verathon Medical (Europe) B.V.
Willem Fenengastraat 13
1096 BL Amsterdam
The Netherlands
Tel: +31 (0) 20 210 30 91
Fax: +31 (0) 20 210 30 92

Verathon Medical (Australia) Pty Limited
Unit 9, 39 Herbert Street
St Leonards NSW 2065
Australia
Tel: 1800 613 603 (Within Australia)
Fax: 1800 657 970
Tel: +61 2 9431 2000 (International)
Fax: +61 2 9475 1201