



# GlideScope Performance. Anywhere.

The GlideScope® Go™ 2 portable video laryngoscope is purpose-built for Emergency Medical Service (EMS) professionals, delivering exceptional visualization, EMS-certified performance, and a design that's intuitive and easy-to-use.



**GlideScope® Go™ 2**  
verathon

# Built for the EMS Environment

## Exceptional Visualization

- Clear visualization on a scratch-resistant 3.5" Corning® Gorilla® Glass<sup>1</sup> screen coupled with anti-fog capabilities
- Optimized image quality and brightness at the vocal cords with Dynamic Light Control (DLC)

## EMS-Certified Performance

- Temperature and vibration-tested for road and air ambulance environments<sup>3</sup>
- Trusted GlideScope® Spectrum™ clinical performance brought to the field

## Intuitive and Easy-to-Use Design

- Seamless workflow and secure connection with magnetic QuickConnect™ technology
- Quality assurance (QA) and quality improvement (QI) made easy with touchscreen controls for auto-recording and media playback<sup>4</sup>

### Rugged and Durable

Corning® Gorilla® Glass<sup>1</sup> and IP67<sup>2</sup> rating for demanding conditions

### 3.5" Touchscreen

Color display on a landscape screen



### Onboard Storage

Up to 40 hours of video storage

### Media Library

Support QA & QI review with snapshots and videos

### Battery Indicator

Quickly see remaining use time

### Rechargeable Battery

Integrated battery provides up to 100 minutes of continual use



Discover the Power of  
GlideScope Go 2 Today

1. Corning and Gorilla are registered trademarks of Corning Incorporated.
2. Ingress protection, IP67, ensures complete dust protection for particles down to 1 mm and submersible up to 1 meter for 30 minutes.
3. IEC 60601-1-12:2014 certified for EMS environments, temperatures of -20 to 50°C and air ambulances including helicopters and fixed wing.
4. Playback available on software versions 1.10 and higher.

verathon

verathon.com

Verathon Inc.  
20001 North Creek Parkway  
Bothell, WA 98011  
USA  
Tel: +1 800 331 2313  
(USA & Canada only)  
Tel: +1 425 867 1348  
Fax: +1 425 883 2896

EC REP

Verathon Medical (Europe) B.V.  
Willem Fenengstraat 13  
1096 BL Amsterdam  
The Netherlands  
Tel: +31 (0) 20 210 30 91  
Fax: +31 (0) 20 210 30 92

Verathon Medical (Australia) Pty Limited  
Unit 9, 39 Herbert Street  
St Leonards NSW 2065  
Australia  
Tel: 1800 613 603 (Within Australia)  
Fax: 1800 657 970  
Tel: +61 2 9431 2000 (International)  
Fax: +61 2 9475 1201